



NEPTUNE 12"

**NEPTUNE HBOT
Hyperbaric Vital Signs Monitor**

- Matchless performance up to a depth of 60 mH2O
- Monitored parameters: ECG 3/5 lead - RESP - NIBP - TEMP (2 channels) - SPO2 - ARR ANALYSIS -ST ANALYSIS
- NIBP technology based on a pneumatic valve supplied with compressed AIR
- Suitable for adults, children and newborns
- 12" tft led colour monitor with high resolution
- Small dimensions and light weight (3.5 kg)
- Special hermetic battery packs with autonomies from 8 hours

**NEPTUNE HBOT
Monitor parametri vitali per camera iperbarica**

- Concepito per funzionare sino a 60 metri di profondità
- Parametri: ECG 3/5 Der - RESPIRO NIBP - 2 CANALI TEMPERATURA - SPO2 - ANALISI ARITMIA - ANALISI ST
- Tecnologia NIBP basata su valvola Pneumatica alimentata ad aria compressa
- Adatto per pazienti adulti, pediatrici e neonati
- Display a led 12" ad alta risoluzione
- Piccole dimensioni a peso contenuto (3,5 kg)
- Gruppo batterie ermetiche con autonomie de 8 ore



**Code n° A56.NNT
NEPTUNE HBOT
MONITOR FOR HYPERBARIC CHAMBER**

FOR FURTHER ACCESSORIES, PLEASE REFER TO THE CURRENT PRICE LIST • PER ALTRI ACCESSORI FARE RIFERIMENTO AL LISTINO PREZZI IN VIGORE